

Colgate® ProClinical® 150 Battery

INSTRUCTIONS

To clean: After use, rinse toothbrush head thoroughly, shake off excess water and store upright.

To replace battery: Firmly hold the bottom of the handle while rotating the top anti-clockwise. Insert 2 new AAA batteries as indicated inside the brush handle and firmly rotate clockwise back into position to close.

To replace brush head: Gently pull the brush head from the handle and press any new Colgate® ProClinical® brush head firmly into place.

Dentists and Hygienists recommend replacing your toothbrush refill every 3 months.

WARNING: Adult supervision required for persons with disabilities. Replace refill if brush head becomes loose due to damage. Battery choke hazard. Do not install batteries backwards, charge, put in fire, or mix with other battery types. May explode or leak causing injury. Replace all batteries at the same time.

Have any questions? Free call 1800 802 307 between 7:00am to 5:00pm Monday to Friday (AEST) or [click here to visit our contact page](#).